

ElectraSpec Installation Files

Please click the link below to

[Download ElectraSpec Installation Files](#)

(Please note: it may take several seconds to start downloading)

After download ESpecInstallation.zip, please COPY all files inside pack to root directory of a USB flash drive (Please refer to Figure 1). Then run installation from flash drive and followed steps given below:

- Go to the flash drive then enter directory \Setups\ESpec;
- Find Setup.exe (Please refer to Figure 2), double click on the file then follow the installation steps
- Please ignore the Subscription option as it is not supported anymore;





Figure 1:

Copy all files and directories to the ROOT directory of a USB drive



If not copy all files and directories to the root directory of a USB drive as show in Figure 3, the installation won't work.

